


OSMO6

Professional reverse osmosis and deionisation unit

HOW IT WORKS

- 1 Pretreatment by activated granular coal** filter 99% of the possible organic matter in suspension, 99% of possible heavy metals, and 99.999% of hypochlorite present in the water.
- 2 The reverse osmosis membrane** rules out 99.999% of organic matter and 94-98% of total dissolved solids.
- 3 Resin filters** remove all the suspended solids up to the trace level.
- 4 Finally the amicrobic filter** ensures good quality Microbiological Water Type II.
The process is complemented by the installation of a UV lamp located inside the storage tank.

The conductivity sensors installed constantly monitor the quality of water produced at every stage.



Production from 240 liters/day

The water produced conforms to the requirements ASTM, CAP, CLSI and ISO 3696 / BS 3997

Don't take any risk with your water

Osmo6 is specifically designed for seminal doses production in boar studs.

The highest quality water for your semen doses

